| | | |
|------------|---|---------|
| D | Gebrauchsanleitung Blutdruckmessgerät..... | (2-13) |
| GB | Instruction for Use Blood pressure monitor | (14-23) |
| F | Mode d'emploi Tensiomètre | (24-34) |
| E | Instrucciones para el uso Tensiómetro | (35-45) |
| I | Istruzioni per l'uso Misuratore di pressione..... | (46-56) |
| TR | Kullanma Talimatı Bilgisayarlı tansiyon ölçer | (57-66) |
| RUS | Инструкция по применению Прибор для измерения артериального давления на предплечье | (67-78) |
| PL | Instrukcja obsługi Ciśnieniomierz | (79-89) |

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Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for the applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

Best regards,
Your Beurer Team

1. Getting to know your instrument

The upper arm blood pressure monitor is used for non-invasive measurement and monitoring of adults' arterial blood pressure.

You can use it to measure your blood pressure quickly and easily, storing the results and displaying the progression of readings together with the average.

A warning is issued for anyone suffering from cardiac arrhythmia.

The values determined are classified and graphically evaluated according to WHO guidelines.

Keep these instructions carefully for further use and also let other users have access to them.

2. Important information

Advice on use

- Always measure your blood pressure at the same time of day, in order to ensure that values are comparable.
- Relax for approx. 5 minutes before each measurement.
- You should wait five minutes before measuring a second time.
- The readings that you take may only be used for information purposes – they are not a substitute for a medical examination! Discuss your results with the doctor. Never use them to make medical decisions independently (e.g. regarding medication and dosage thereof).
- There may be incorrect measurements where there is cardiocirculatory disease, as well as with very low blood pressure, circulatory disorders, dysrhythmia and other preexisting diseases.
- Only use the instrument on persons whose upper arm has the right circumference for the instrument.
- You can operate the blood pressure monitor with both batteries and mains. Note that it is only possible to save data when your blood pressure monitor is receiving power. As soon as the batteries wear out or you disconnect the mains adapter from the mains, the blood-pressure instrument loses the date and time.

- The automatic switch-off function switches off the blood pressure monitor in order to preserve the batteries, if no button is pressed within one minute.

Storage and Care

- The blood pressure monitor is made up of precision electronic components. Accuracy of readings and the instrument's service life depend on careful handling.
 - You should protect the device from impact, moisture, dirt, major temperature fluctuations and direct exposure to the sun's rays.
 - Never drop the device.
 - Do not use near strong electromagnetic fields, i.e. keep it away from any radio systems and mobile phones.
 - Only ever use the cuffs provided with the monitor or original replacement cuffs. Otherwise erroneous results will be recorded.
- Do not press any buttons until the cuff is in position.
- If the instrument is not used for any length of time, we recommend removing the batteries.

Advice on batteries

- Batteries can be fatal if swallowed. You should therefore store the batteries and products where they are inaccessible to small children. If a battery has been swallowed, call a doctor immediately.

- Batteries should not be charged or reactivated with any other means, nor should they be taken apart, thrown in the fire or short-circuited.
- Remove the batteries from the instrument if they are worn out or if you are not going to use the instrument for any length of time. This prevents any damage as a result of leakage. Always replace all the batteries at the same time.
- Never use different types of battery, battery brands or batteries with different capacities. You should preferably use alkaline batteries.


Repair and disposal


- Batteries do not belong in domestic refuse. Used batteries should be disposed of at the collection points provided.
- Never open the instrument. If these instructions are not heeded, the warranty will be null and void.
- Never attempt to repair the instrument or adjust it yourself. We can no longer guarantee perfect functioning if you do.
- Repairs may only be performed by Beurer Customer Service or authorized dealers. However, always check the batteries and replace them if necessary prior to making any complaint.
- The appliance should be disposed of according to Regulation 2002/96/EC-WEEE (Waste Electrical and Electronic Equipment). In case of queries, please contact the municipal authorities responsible for waste disposal in your area.

3. Unit description

1. Cuff tube
2. Cuff
3. Cuff connector
4. Touch screen activation switch
5. USB interface
6. Cuff holder
7. Cuff connector port (left side)
8. Memory button MEM
9. START/STOP button
10. Mains adapter port
11. Display

Touch screen activation switch

The device features a touch screen display. To prevent the screen from being activated accidentally, keep the touch screen activation switch set to "OFF" when the device is not in use. To operate the device, move the touch screen activation switch to "ON". When the touch screen display is touched ("START/STOP" button or "MEM" button), a beep is output.


Tip: You can switch off the device at any time by setting the touch screen activation switch to the "OFF" position.

Cuff holder

Thanks to the cuff holder on the underside of the device, you have the option of conveniently storing the cuff. To do so, push out the cuff holder with both thumbs until it locks into place (see also Chapter: Inserting battery).

Icons in the display:

1. Time and date
2. Systolic pressure
3. Diastolic pressure
4. Measured pulse
5. Pump up, release air (arrow)
6. Memory display, average value (R), morning (R_M), evening (P_M)
7. WHO classification

8. „Change battery“ icon 
 9. Arrhythmia recognition 


PC kit

With the additionally available accessory set 655.17, consisting of a CD with PC software and a data-transfer cable, you can evaluate your data extremely effectively. You will find the installation instructions in the PC kit.

4. Prepare measurement

Inserting battery

- Push the cuff holder upwards at the rear of the device until it locks into place.
- Open the battery compartment lid
- Insert four AA 1.5 V alkaline batteries. making absolutely sure that you insert them with the correct polarity as marked. Never use rechargeable batteries.
- Replace the battery cover carefully.

If the battery change „“ is continuously illuminated, measurement is no longer possible and you must replace all

the batteries. Once batteries have been removed from the device, the time must be reset.

Used batteries do not belong in the household waste. You are legally obliged to dispose of the batteries. Dispose of them via your specialist electrical supplier or local collecting point for recyclable waste.

Note: Batteries containing pollutant substances are marked as follows: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

Setting date and time

You should set the date and time without fail. Otherwise, you will not be able to save your measured values correctly with a date and time and to access them again later.

The time is displayed in the 24-hour format.

Tip: If you press and hold the "MEM" button, you can set the values more quickly.

To set the date and time, proceed as follows:

- Move the touch screen activation switch to the "ON" position.
- Press the "START/STOP" and "MEM" buttons simultaneously. The year display will start to flash. Set the year with the "MEM" button and confirm with "START/STOP".
- Then set the month, day, hour and minutes and confirm each setting with "START/STOP".

- Pressing the "START/STOP" button again will switch the display off.

Mains operation


You can also operate this device with a mains adapter. No batteries should be in the battery compartment for this. The mains adapter is available from retailers or from the service address under order no. 071.19. The blood pressure monitor should only be used with the mains adapters specified here. The mains adapter should only be connected to the mains voltage indicated on the rating plate. Once you have disconnected the device from the mains, the date and time on the blood pressure monitor are lost. Any measurement results that have been stored are however retained.

5. Measuring blood pressure

Please ensure the device is at room temperature before measuring.




Positioning cuff

Fit the cuff round your bare left upper arm. Blood circulation in the arm should not be restricted by tight clothing or other objects. The cuff should be placed on the upper arm so that the lower edge is 2 to 3 cm above the bend of the




elbow and above the artery. The tube should be in line with the centre of the palm.

Now place the free end of the cuff snugly, but not too tightly, around the arm, and fix it with the Velcro fastener. The cuff should be fitted tight enough to allow just two fingers to fit beneath the cuff. Insert the cuff tubing into the socket for the cuff attachment.

Important: The instrument should only be operated with the original cuff. The cuff is suitable for an arm circumference of 22 to 30 cm. A larger cuff for upper arm circumferences from 30 to 42 cm is available from retailers or the service address under order no. 162.796.

Correct posture

- Rest for approx. 5 minutes before each measurement. Otherwise there may be divergences.
- You can perform the measurement either sitting or lying down. Always make sure that the cuff is on a level with your heart.
- In order not to distort the result, it is important to keep still during the measurement and not talk.

Measuring blood pressure

- Move the touch screen activation switch to the "ON" position.
- Put on the cuff as described previously and assume the position in which you want to carry out the measurement.
- Switch on the blood pressure monitor with the "START/STOP" button. After the full-screen display, the most recently used user memory appears ( or ). To change the user memory, press the "MEM" button and confirm your selection with the "START/STOP" button. If no button is pressed, the most recently used user memory is automatically used after 5 seconds.
- Before the measurement, the last saved test result is briefly displayed. If there is no measurement in the memory, the instrument always displays the value "0".
- The cuff is pumped up to 190 mmHg. Cuff air pressure is released slowly.
If a tendency towards high blood pressure is already detectable, the cuff is pumped up again and cuff pressure

increased further. As soon as a heart rate is detected, the heart rate symbol  is displayed.

- Heart rate, systolic and diastolic blood pressure are displayed.
- You can cancel the measurement at any time by pressing the "START/STOP" button or by pushing the touch screen activation button to "OFF".
- "Er" appears if it has not been possible to perform the measurement properly. Observe the section in these instructions on error messages/troubleshooting and repeat the measurement.
- The test result is saved automatically.
- To switch off the device, press the "START/STOP" button or push the touch screen activation switch to "OFF". If you forget to switch off the device, it switches off automatically after approx. 1 minute.

Wait at least 5 minutes before taking another measurement!

6. Evaluating results

Cardiac arrhythmia:

This instrument can identify possible cardiac arrhythmia disorders during measurement and if necessary indicates the measurement with the flashing icon .

This may be an indicator for arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heart beat. The symptoms (omitted or premature heart beats,

slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor.

Repeat the measurement if the flashing icon  is displayed after the measurement. Please note that you should rest for 5 minutes between measurements and not talk or move during the measurement. If the icon  appears often, please contact your doctor. Any self-diagnosis and

treatment based on the test results may be dangerous. It is vital to follow your doctor's instructions.

WHO classification:

The World Health Organization (WHO) and National High Blood Pressure Education Program Coordinating Committee have developed a blood pressure standard, according to which areas of low- and high-risk blood pressure are identified. This standard, however, is a general guideline as individual's blood pressure varies among different people and different age groups...etc.

| Range of blood pressure values | Systolic (in mmHg) | Diastolic (in mmHg) | Measure |
|--------------------------------|-----------------------|------------------------|-------------------------------------|
| Grade 3: Severe hypertension | >= 180 | >= 110 | Seek medical advice |
| Grade 2: Moderate hypertension | 160–179 | 100–109 | Seek medical advice |
| Grade 1: Mild hypertension | 140–159 | 90–99 | Have it checked regularly by doctor |
| High-normal | 130–139 | 85–89 | Have it checked regularly by doctor |
| Normal | 120–129 | 80–84 | Check it yourself |
| Optimal | <120 | <80 | Check it yourself |

Source: WHO, 1999

It is important that you consult with your physician regularly. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk.

The WHO classification in the display shows the range of the measured blood pressure. If the values for systolic and diastolic pressure are in two different WHO ranges (e.g. systolic in the grade hyperten-

sion range and diastolic pressure in the normal range), the WHO classification on the unit always indicates the higher range (grade hypertension in the example described).

7. Saving, retrieving and deleting results

- The results of each successful measurement are stored together with date and time. With more than 60 items of measured data, the earliest items of data measured are lost.
- Move the touch screen activation switch to the "ON" position.
- Make a selection with the "MEM" button, and then confirm the desired user memory with the "START/STOP" button. If you press the "MEM" button again, the average value " \bar{P} " of all the stored measured values in the user memory will be displayed. If you press the "MEM" button again, the average value of the last 7 days for the morning measurement will be displayed (Morning: 5 a.m. - 9 a.m., display " \bar{P}_m "). If you press the "MEM" button again, the average value of the last 7 days for the evening measurement will be displayed (Evening: 6 p.m. - 8 p.m., display " \bar{P}_e "). If you continue to press the "MEM" button, the most recent individual measured values with date and time are displayed in turn.
- You can delete the memory by pressing and holding the "MEM" button for 3 seconds. All the values in the current user memory are deleted after three beeps are output.

- To switch off the device, press the "MEM" button again or the "START/STOP" button or push the touch screen activation switch to "OFF".
- If you forget to switch off the device, it will switch off automatically after 2 minutes.

8. Error messages/trouble-shooting

In case of faults, the "Er_" message appears in the display. Error messages can occur when

- the blood-pressure values are unusually high or low (Er1 or Er2 appears in the display),
- the cuff is too tight or too loose (Er3 or Er4 appears in the display),
- the inflation pressure is higher than 300 mmHg (Er5 appears in the display),
- inflation takes longer than 3 minutes (Er6 appears in the display),
- there is a system or device error (ErA, Er0, Er7 or Er8 appears in the display),
- the batteries are almost empty .

In such cases, repeat the measurement. Ensure you do not move or speak.

If necessary, re-insert or replace the batteries.

9. Cleaning and storing the instrument

- Clean your blood pressure computer carefully only with a slightly moistened cloth.
- Do not use detergents or solvents.
- On no account must you immerse the computer in water, otherwise liquid can enter it and cause damage.
- When storing the device, make sure that no heavy objects are placed on top of it. Remove the batteries. The cuff tube should not have any sharp kinks.

10. Specifications

| | |
|----------------------------------|---|
| Model No. | BM 58 |
| Method of Measurement | Oscillometric, non-invasive blood pressure measurement on the upper arm |
| Range of Measurement | Pressure 0–300 mmHg, systolic 30–280 mmHg, diastolic 30–280 mmHg, pulse 40–180 beats/minute |
| Accuracy of the pressure reading | systolic ± 3 mmHg / diastolic ± 3 mmHg / pulse $\pm 5\%$ of the displayed value |
| Measuring uncertainty | max. permissible standard deviation according to clinical testing: systolic: 8 mmHg / diastolic: 8 mmHg |
| Sets of Memory | 2 x 60 memory places |

| | |
|-----------------------|--|
| Unit Dimension | L100 x W150 x H58 mm |
| Unit Weight | Approx. 364 g |
| Cuff Size | 22 – 30 cm |
| Operating Temperature | +5 °C to +40 °C, < 90 % RH |
| Storage-Temperature | -20 °C to +55 °C, < 95 % RH |
| Power Supply | 4 x 1,5 V AA batteries (alkaline type LR6) |
| Battery Life | For approx. 500 measurements, depending on blood pressure level or inflation pressure |
| Accessories | With storage pouch, instruction manual, 4 pcs “AA” battery |
| Classification | Application part Type B |
| Protection class | Internal power supply, IPX0, no AP or APG, continuous operation |
| Key to symbols: | Application part Type B  |
| | Important! Read the operating instruction!  |
| | Manufacturer  |

These specification are subject to change without notice for purpose of improvement.

11. Adapter

| | |
|------------|---|
| Model No. | FW 7575M/EU/6/06 |
| Input | 100–240 V, 50–60 Hz |
| Output | 6 V DC, 600mA, only in connection with beurer blood pressure monitor |
| Supplier | Friwo Gerätebau GmbH |
| Protection | This device is double insulated and protected against short circuit and overload by a primary thermal fuse. Make sure to take the batteries out of the compartment before using the adapter |

double insulated / equipment class 2

| | |
|----------------------------------|---|
| Enclosures and Protective Covers | Equipment enclosed to protect against contact with live parts, and with parts which can become live (finger, pin, hook test). The operator shall not contact the patient and the output plug of AC adaptor simultaneously. |
|----------------------------------|---|

- This unit is in line with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility (EMC). Please note that portable and mobile HF communication systems may

interfere with this instrument. For more details, please contact customer service at the address indicated.

- The instrument is in line with the EU Medical Products Directive 93/42/EC, the German medical products act and European standards EN1060-1 (non-invasive blood pressure monitors, Part 1: general requirements) and EN1060-3 (non-invasive blood pressure monitors, Part 3: Supplementary Requirements for Electromechanical Blood Pressure Measurement Systems).

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Producent Friwo Gerätebau GmbH