

Productspecification

Art.-No.: REF 63101

Article: LEINATEX – First Aid Dressing 60 x 80 cm

<u>Item:</u> Surgical sheet with aluminium coating individually in blister

Packing sterilized, white, without optical brightener

Thickness: 1mm

Breaking strength: 120 N / 5 cm

Material: absorbent wad: non woven, knit type

Raw material: viscose rayon 70%

Polypropylene fibbers 30%

Weight per unit area: 150g/m² + 10g/m²

Topside cover: net shaped polyethylene film metallised

with aluminium $15 \text{ g/m}^2 \pm 5 \text{g/m}^2$

Backside cover: PP-non woven 25g/m² ± 5g/m²

(both top-side and back-side hydrophob pervious to

secretion of wound)

Size: BR: $(600 \pm 15) \times (800 \pm 20) \text{ mm}$

Folding size: $100 \times 155 \times 20 \pm 5 \text{ mm}$

Water absorption capacity: minimum 800g/m² (test method according to DIN 53923)

This item is suitable for burned wound and does not stick to the wound.