## ri-focus® LED

## WIRELESS HIGH PERFORMANCE LED HEADLIGHT POWERED BY LI-BATTERIES OR RECHARGEABLE BATTERIES

## green efficiency

- ri-focus® LED With 1 set of CR 123 A lithium batteries Art. No.: 6090
- With 1 set of NiMH rechargeable AAA batteries and plug-in charger

230 V Art. No.: 6091 120 V Art. No.: 6092

Plug in chargers 230 V Art. No.: 11284 120 V Art. No.: 11286

Batteries Lithium CR 123 A 2 pcs Art. No.: 11287 AAA NiMH rechargeable batteries 4 pcs Art. No.: 11288

Foam-rubber band 120 mm long Art. No.: 11289

Fabric points Ø 35 mm Art. No.: 11293 Ø 47 mm Art. No.: 11294

Battery compartment cover For lithium CR 123 A Art. No.: 12680 For AAA NiMH rechargeable batteries Art. No.: 12681

- One hand operation
- Wireless and balanced headband with built in battery compartment
- More efficient diagnosis due to white high performance LED (140 lumen)
- Service life of LED up to 50.000 hours
- Significantly less energy consumption and heat generation
- Headband version, with built-in battery compartment and integral charge jack, can be optionally fitted with 2 lithium batteries size CR 123 A, 4 batteries size AAA (alkaline) or 4 rechargeable batteries size AAA (NiMH).
- Headband cleaning made easy by inner, removable and washable padding.
- Particularly comfortable with balanced, infinitely adjustable headband.
- On/off switch on battery compartment.
- Charge jack for plug-in charger integrated in battery compartment.
- Focusable LED lamp at a distance of 40 cm /15.75 inch approx. Ø 40 mm /1.57 inch Ø 200 mm / 7.87 inch, infinitely adjustable.
- White LED with 6.500 Kelvin.
- Battery operation: approx. 4 hrs. with lithium batteries, approx. 90 min. with alkaline batteries Type AAA.
- Rechargeable battery operation: approx. 90 min. with fully charged battery.
- Weight: 210 g / 0.46 lb without batteries, 242 g / 0.54 lb with CR 123 A lithium batteries, 259 g / 0.57 lb with rechargeable AAA batteries.



| Distance               | 400 mm / 15.75 inch | 200 mm / 7.87 inch |
|------------------------|---------------------|--------------------|
| Lux max                | 10.000              | 23.000             |
| Circle of light max. Ø | 200 mm / 7.87 inch  | 120 mm / 4.73 inch |
| Circle of light min. Ø | 40 mm / 1.57 inch   | 30 mm / 1.18 inch  |

## ri-focus® LED