

BLOC CONNECTOR

Description

Accessory device with a cone-shaped tip at one end and a Luer Lock connector at the other end.

Class I sterile, medical device, complying with Directive 93/42/EEC.
Single use.
Latex-free.
Individual packaging.
Sterilised by ethylene oxide.

Indication

The Bloc Connector is used to connect any device fitted with a Luer tip (syringe, needle, catheter...) to a central vacuum system or to suction tubes.

Instructions for use

- Connect the cone-shaped tip to the suction tube.
- Connect the cylinder section to the Luer tip.
- Rotate the cylinder section to lock it onto the Luer tip and provide a perfect seal between the two.

Reference

1990300 BLOC CONNECTOR *box of 25 units.*

Storage

To store in a dark place.