



Little Junior QCPR

Providing feedback based on just watching students perform can be challenging for instructors. The next generation of the Little Junior manikin includes QCPR measurement and feedback technology that takes the guesswork out of instruction all from the palm of your hand. Little Junior QCPR works with an Instructor app that allows you to monitor and coach up to six students at once, and a Learner app that lets students monitor their own performance.

Little Junior QCPR helps instructors improve CPR training quality and efficiency in their courses, while adding a fun and engaging element for students that promotes learning and motivation. (1)

All existing Little Junior manikins can be quickly and easily upgraded to QCPR with the Little Junior QCPR Upgrade Kit.





Little Junior QCPR

Feedback Options

Little Junior QCPR's feedback technology provides:

- Real-time feedback on compressions and ventilations
- A performance summary with an overall score and improvement suggestion
- · Details on compression release, depth and rate, ventilation volume, and number of compressions/ventilations/cycles

Training based on the use of QCPR technology guided by an instructor has been shown to be superior to instructor-based feedback training for chest compression technical skill acquisition.(2)







QCPR Instructor App

For instructors to monitor multiple students

- Free download on App Store
- Connects to 1-6 manikin(s)
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, postsession
- Results are saved in the app
- Race mode



QCPR Learner App

For students to monitor their own performance

- · Connects to I manikin
- Shows live and summative feedback of CPR performance
- Provides suggestions for improvement, post-
- Results are saved in the app
- Free download on App store and Google Play





SkillGuide

- Provides a tethered, "plug and play" option
- Can be used standalone or as a supplement to either of the apps

ORDERING INFORMATION

128-01050 Little Junior QCPR 128-03050 Little Junior QCPR Dark 129-01050 Little Junior QCPR 4-Pack 129-03050 Little Junior QCPR 4-Pack Dark 128-60750 Little Junior QCPR Upgrade Kit

Accessories and Consumables

123-30050 SkillGuide with extension cable 123-50550 Extension Cable for SkillGuide 183210 Junior Airways (pkg. 25) 123211 Junior Airways (pkg. 100) 15120103 Manikin Face Shields (6 rolls of 36 ea) Manikin Wipes (pkg. 50) 152400 152401 Manikin Wipes (pkg. 1200) 183010 Junior Faces (pkg. 6) 183015

Junior Faces Dark (pkg. 6)



Visit Laerdal.com/Little|uniorQCPR to learn more about Little |unior QCPR and how you can upgrade your existing Little Junior manikins.

©2018 Laerdal Medical. All rights reserved. 18-16720 Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google LLC. 1. Smart et al. Does real-time objective feedback and competition improve performance and quality in manikin CPR training - a prospective observational study from several European EMS. Scandinavian Journal of Trauma, Resuscitation and Emergency Medicine. 2015 2. Cortegiani et al. Use of a Real-Time Training Software (Laerdal QCPR®) Compared to Instructor-Based Feedback for High-Quality Chest Compressions Acquisition in Secondary School Students: A Randomized Trial. PLoS One. 2017

